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Ziegelei 1
72336 Balingen-Frommern
Germany

www.kern-sohn.com

+0049-[0]7433-9933-0

+0049-[0]7433-9933-149

info@kern-sohn.com

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The screenshot shows the USR IOT web interface. The top navigation bar includes the USR IOT logo and the text 'Communication Expert of Industrial IoT'. A sidebar on the left lists menu items: Status, Network, Port, Gateway, Cloud Service, and System. The main content area is divided into three sections: 'Status', 'Network', and 'Port'. The 'Status' section contains a table with system information. The 'Network' section contains a table with network configuration details. The 'Port' section contains a table with port status and traffic statistics.

System	
Model Name	USR-N520
Firmware Version	V2.0.10
Type	H7-6
Run Time	0day: 5hour: 50min

Network	
MAC Address	D4-AD-20-AE-C7-40
Current IP Address	192.168.0.7
Preferred DNS Server	114.114.114.114
Alternate DNS Server	223.5.5.5

Port	
Status of Port	Port1
Conn Status A(ETH)	CONNECTED(1)
TX Count A(ETH)	0 bytes
RX Count A(ETH)	8409 bytes
Conn Status B(ETH)	IDLE
TX Count B(ETH)	0 bytes

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The screenshot shows the USR IOT web interface. On the left, a navigation menu is visible with the following items: Status, Network (expanded), IP Config (selected), Port, Gateway, Cloud Service, and System. The main content area is titled 'IP Configuration' and 'IP Configuration of WAN Port'. Under the 'Configuration' section, the following settings are displayed:

Method of IP Obtaining	Static IP
DNS	Auto
IP Address	192.168.0.7
Subnet mask	255.255.255.0
Gateway	192.168.0.1
Preferred DNS Server	114.114.114.114
Alternate DNS Server	223.5.5.5

A 'Save&Apply' button is located at the bottom right of the configuration area.

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5.6 Serielle RS232 Schnittstellen

The screenshot shows the 'USR IOT' web interface. On the left is a navigation menu with 'Port' selected. The main content area is titled 'UART TO ETH' and 'Data transmission parameter configuration'. Under the 'SETTING' header, the 'Port' tab is active. The configuration includes: Baud rate (115200), Data bits (8), Parity (None), Stop bits (1), Serial Mode (RS485), Current Serial Mode (RS485), Flow ctrl (NONE), UART Packet Length (0), UART Packet Time (0), Sync Baudrate(RFC2217) (ON), and an unchecked 'Enable Uart Heartbeat' checkbox. A 'Save&Apply' button is at the bottom right.

The screenshot shows the 'USR IOT' web interface with the 'Socket' tab active. The main content area is titled 'SOCKET A' and 'Data transmission parameter configuration'. Under the 'SETTING' header, the 'Socket' tab is active. The configuration includes: Working Mode (TCP Server), Maximum Sockets supported (8), Local Port Number (23), PRINT (OFF), Modbus Poll (unchecked), and an unchecked 'Enable Net Heartbeat' checkbox. Below this is the 'SOCKET B' section with 'Operating Mode' set to 'None'. A 'Save&Apply' button is at the bottom right.

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