

# Digital Refractometer



# **Operating Manual**

1.	General information	3
2.	Introduction	4
3.	Display & operating buttons	5
4.	Preparing before operating	6
5.	Turn on & measure	6
6.	Calibration	8
7.	Changing scale & temperature unit	1(
8.	Turning off	11
9.	Cleaning & maintenance	1
10	Disposal	11
11	· Technical data	12
12	. Error codes	12
13	. Models and scales	13
<u>/</u> !	7	

Conten

Carefully read through the operating manual even if you have prior experience with KERN refractometers.

### 1.1 Intended use

The refractometer is a measuring instrument for determining the refractive index of transparent substances in liquid or in some cases also in the solid state. It is used to observe the behaviour of light as it passes from a prism with known properties to the substance being tested. Use of the refractometer for other purposes is contrary to its intended use and may be hazardous. The manufacturer shall not be liable for any damages caused by improper use.

## 1.2 Warranty

The warranty shall be void in the event of Failure to observe the instructions in the operating manual Use for purposes other than those described Modifications or opening the device housing Mechanical damage and/or damage resulting from media, liquids, natural wear and tear

This digital refractometer cannot measure any liquid that is highly corrosive to metal or glass. When measuring liquids that are corrosive to plastics or react chemically with plastics, be careful not to drop the measured liquid onto the shell. Otherwise it will corrode the shell.





### Note ·

1. When used outdoors, please avoid strong iaht so

as not to affect the measurement accuracy 2 Please keep the instrument in a stable and still statement and position.

After turning on, clean the sample tank with distilled water and then dry it. Now fill the sample up to the mark, close the cover and press "READ".

Measurement HHH

measurement, press

"READ" for 2 seconds

.....

204T

LLL





5.3 Average value measurement

Press "READ" for 2 seconds. The device starts an automatic measurement series of 15 measurements and shows the average value. Afterwards, the device automatically turns back to the normal measuring mode.



Remaining measurements

The refractometer can only be calibrated with distilled water. To do this, fill the sample tank with distilled water up to the mark and close the cover.



While "CAL" is flashing in the display, press "CAL" again to start the calibration. When the calibration is finished, the display shows "End". After approx 10 seconds the device automatically returns to normal mode.



If the calibration was not completed successfully. an error code appears in the display. Here, for example, A01.



Further error codes can be found in the appendix.

We recommend calibrating the refractometer.

- when commissioning
- after a strong shock
- after longer transport
- after a change of location with a large temperature difference
- If the device has not been used for a long time

Always use distilled water and make sure that the refractometer, the water and the environment are at the same temperature.

## Changing scale & termperature unit

## 7.1 Changing scale



## 7.2 Changing termperature unit



If exceed the temperature limitations, the signs "HHH" or "LLL" would show.





If without any operations for 1 minute. the instrument would be automatically turned off.

- 1. To avoid damages to the prism and the sample tank, clean them with distilled water after each
- 2 Drv it with a soft cloth afterwards
- 3. Do not use hard or abrasive objects for cleaning.
- 4. Do not leave any residue in the sample tank.
- 5. If the refractometer is not going to be used for a longer time, remove the battery and store it at a cool and dry place.

The packaging consists of environmentally friendly materials which can be disposed of via local recycling facilities. The device and storage box should be disposed

code	In
A01	Beyond the temperature
A02	During calib or solution v
A03	This instrum failure.

1. Technical data

Scale + accuracy + resolution

Temperature

Automatic Temperature Compensation

Lifetime of the battery

Net weight

Minimum sample volume

Averaging measurement

verall dimensions L×W×H

## of by the operator in accordance with applicable national or regional regulations at the place of

NOTE: In accordance with the Battery Ordinance (BattV) batteries must not be disposed of in householdwaste. The end user is legally obliged to return them.

1	Depents to the model
	0,0-40,0 °C/32,0-104,0 °F
	Yes
	0.2 - 0.3 ml (Marking ring)
	60 seconds
	15 measurements
	1 × AAA 1.5 V
	Approx. 10.000 measurements
I	125×65×30 mm
	140 g (without battery)

## structions

e scope of calibration re. (0.0℃~40.0℃) 

libration, no solution wrona. 5

ment has a hardware

## 3. Models and scales

-	Model	Scale	No.	Range	Linit	Resolution	Accuracy
-	ORM 50BM	Brix	S01	0.0~50.0	%	0.1%	±0.2%
	ORIM SOBM	Refractive Index	S02	1.3330~1.4200	00	0.0001nD	+0.00030D
	ORM 1RS	Brix	801	0.0-90.0	%	0.1%	±0.2%
	0.000 1100	Refractive Index	S02	1.330-1.5177	nD	0.0001nD	±0.0003nD
-	ORM 1SU	Fructose	S01	0.0~68.9	%	0.1%	±0.2%
	ORM 1SU	Glucose	502	0.0-59.9	%	0.1%	±0.2%
		Brix	303	0.0-90.0	%	0.1%	±0.2%
8		Refractive Index	804	1.3330~1.5177	nD	0.0001nD	±0.0003nD
ructose	ORM 25U	Lactose	801	0.0~16.5	%	0.1%	±0.2%
£.	010/1250	Mallose	S02	0.0~15.6	%	0.1%	±0.2%
		Dextran	503	0.0~10.6	%	0.1%	±0.2%
		Brix	804	0.0-50.0	~	0.1%	±0.2%
-	ORM 1HO	Honey Water	801	5.0-38.0	%	0.1%	±0.2%
5	ORM THO	Honey Baume	802	33.0-48.0	*Bé	0.1	+0.2
Honey		Brix	503	0.0-90.0	%	0.1%	±0.2%
-		Refractive Index	804	1.3330-1.5177	nD	0.0001nD	±0.0003nD
-	ORM 1NA	Salinity (NaCl) %	801	0.0-28.0	%	0.1%	+0.2%
	ORM INA	Salinity (NaCl) %	802	0-280	~ ~	1%	+2%
		Specific Weight	S03	1.000-1.220		0.001	±0.002
		Brix	S04	0.0-50.0	%	0.1%	±0.2%
£.		Refractive Index	805	1.3330~1.4200	nD	0.0001nD	±0.0003nD
Fee Sec	ORM 1SW	Salinity Seawater	S01	0-100	Se	1%	±0.0003nD ±2%
ಹ	Orea 18W	Chlorinity Seawater	S02	0~57	744 564	1%e	±2%
		Specific Weight	502	1.000~1.070		0.001	±0.002
		Brix	S04	0.0~50.0	· %	0.001	±0.002
		Refractive Index	805	1.3330~1.4200	nD	0.0001nD	±0.0003nD
-	ORM 1AL	Alcohol Mass.	801	0-72	11D	1%	±1%
5	ORM 1AL	Alcohol Vol.	S02	0-80	<u>n</u>	1%	±1%
4 ochol		Brix	502	0.0-50.0	%	0.1%	±0.2%
ž		Refractive Index	505	1.3330~1.4200	nD	0.0001nD	±0.0003nD
-	ORM 1BR	Plato	801	0.0-30.5	*P	0.000 mD	±0.3
2	Orea IBR	SG Wort	802	1.000-1.130		0.001	±0.002
Beer		Brix	503	0.0-50.0	~	0.1%	±0.2%
~		Refractive Index	804	1.3330~1.4200	nD	0.0001nD	±0.0003nD
-	ORM 1WN	Oechsle	S01	0-150	*Oe	1	±2
	OPMI IVIN	Vol%	802	0.0-22.0	- Cu	0.1%	+0.2%
		KMW (Babo)	502	0.0-25.0		0.1 %	±0.2
•		Brix	804	0.0-50.0	~	0.1%	±0.2%
/wo	ORM 2WN	Oechsle France	S01	0~230	"Oe	1	+2
-	ORM 200N	Vol%	502	0.0~22.0	- Cu	0.1%	±0.2%
		KMW (Babo)	502	0.0-25.0		0.1	±0.2
		Brix	S04	0.0-50.0	%	0.1%	±0.2%
-	ORM 1CO	Coffee TDS 1	801	0.0~25.0		0.1 %	±0.2
	ORM 100	Brix	502	0.0~50.0	%	0.1%	±0.2%
0		Refractive Index	502	1.3330~1.4200	nD	0.0001nD	±0.0003nD
Coffee		Coffee TDS 2	801	0.00-25.00	112	0.001	±0.20
õ	ORM 2CO	Brix	S02	0.00-30.00	5	0.01%	+0.20%
		Refractive Index	802	1.3330~1.4200	nD	0.0001nD	±0.0003nD
-	-	Urine Human	801	1.000-1.050	110	0.001	±0.002
	ORM 1UN	Serum Protein	502	0.0~12.0	g/100ml	0.001	±0.2
		Brix	502	0.0-50.0	S.	0.1%	±0.2%
		Refractive Index	804	1.3330-1.4200	nD	0.0001nD	±0.0003nD
er H		Urine Cat	S04	1.000-1.060		0.0001110	±0.000311D
~	ORM 2UN	Urine Dog	802	1.000-1.060		0.001	+0.002
		Brix	S02	0.0~50.0	· %	0.001	±0.002
		Refractive Index	804	1.3330~1.4200	nD	0.0001nD	±0.0003nD
-		Cleaner	S01	(-60.0)=0.0	10 10	0.1°C	±0.5°C
	ORM 1CA	AdBlue®	S02	0.0-51.0	%	0.1%	±0.5 C
		Battery Fluid	S02	1.000~1.500	-	0.001	±0.2% ±0.005
£.		Battery Fluid Brix	S04	0.0~50.0	· %	0.001	±0.005
		Brix Refractive Index	804	1.3330-1.4200	% nD	0.1% 0.0001eD	±0.2% ±0.0003nD
2		Ethylenglycol (%)	S05	1.3330-1.4200	10 %	0.0001nD	±0.0003nD ±0.5%
Indu					70 10	0.1% 0.1°C	±0.5%
r / Indu							
Car / Industry		Ethylenglycol (*C)	S02	(-50.0)=0.0			
Car / Indu	ORM 2CA	Ethylenglycol (*C) Propylenglycol (%) Propylenglycol (*C)	S02 S03	(-50.0)=0.0 0.0=100.0 (-60.0)=0.0	% *C	0.1% 0.1%	±0.5% ±0.5°C

2024/09 V3.2