

PRODUCT BROCHURE

KERN 336-37

Date: 19.03.2025

ENGLISH



Individual weight, OIML Class F2, 10 kg, Block, stainless steel

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

- Individual weight, OIML Class F2, 10 kg, Block, stainless steel
- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section
- For the verification you need a box. Please add a box for single weights.

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	336-37
Model-Serie	336
Assortment	KERN
Product Type	Test weights
Product group	OIML F2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761232815
Customs tariff number	84239090

Technical data - Primary

Nominal value	10 kg
Tolerance (OIML)	160 mg
Case aluminium	346-070-600
Case plastic	not available
Case wood	not available

Technical data - Construction

Dimensions	193×96×117 mm
Material	stainless steel

Technical data - Measuring system

Nominal value	10 kg
OIML class	F2
OIML conform	yes
Tolerance (OIML)	160 mg

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	193×96×117 mm
Gross weight	10 kg
Net weight	10 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-444
-------------------	---------

FUNCTIONS

Standard



ACCESSORIES

Model	Description
346-070-600	Aluminium Protective Box KERN 346-070-600

SERVICES

Model	Description
962-444	DAkkS Calibration KERN 962-444
962-444R	DAkkS Recalibration KERN 962-444R
920-30	Marking with stroke number KERN 920-30
921-444	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-444
920-301	Marking with stroke number (Last digits of serial number) KERN 920-301
920-321	Marking by etching (last digits of the serial number) KERN 920-321

RELATED ARTICLES DATA

Model name	Nominal value	Tolerance (OIML)	Case aluminium	Case plastic	Case wood
336-36	5 kg	80 mg	346-060-600	not available	not available
336-37	10 kg	160 mg	346-070-600	not available	not available
336-38	20 kg	300 mg	346-080-600	not available	not available
336-39	50 kg	800 mg	346-090-600	not available	not available