

PRODUCT BROCHURE

KERN 303-056

Date: 10.03.2026

ENGLISH



Weight set, E1, 1 mg - 500 g, Knob, stainless steel polished, in aluminium case


CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information



ENGLISH

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

Model article-code	303-056
Model-Serie	303
Assortment	KERN
Product Type	Test weights
Product group	OIML E1 Test weights
Product usage type	Main article
EAN code (Model)	4045761223035
Customs tariff number	84239090
Country of Origin	CN

Nominal value	1 mg - 500 g
---------------	--------------

Case aluminium	313-050-600
Weight set case	in aluminium case

Material	stainless steel polished
Is weight set	yes
Type of case	in aluminium case

Nominal value	1 mg - 500 g
OIML class	E1
OIML conform	yes

Dimensions packaging (W×D×H)	210×140×98 mm
Gross weight	2 kg
Net weight	1.6 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	400 g
Delivery time	1 d

Optional Services

DAkkS Calibration	963-204
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
<u>963-204</u>	DAkkS Calibration KERN 963-204
<u>962-204R</u>	DAkkS Recalibration KERN 962-204R
<u>920-713</u>	Marking (Last digits of serial number) KERN 920-713

SPARE PARTS

Model	Description
<u>313-050-600</u>	Aluminium Protective Box KERN 313-050-600
<u>315-242</u>	Tweezers KERN 315-242
<u>318-270</u>	Dusting Brush KERN 318-270