

PRODUCT BROCHURE

KERN 334-04

Date: 09.03.2026

ENGLISH



Weight set, F2, 1 g - 200 g, Knob, stainless steel fine turned, in wooden case

CONTACT


KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	334-04
Model-Serie	334
Assortment	KERN
Product Type	Test weights
Product group	OIML F2 Test weights
Product usage type	Main article
EAN code (Model)	404576 10059 14
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 g - 200 g
Case wood	353-440-200
Weight set case	in wooden case

Technical data - Construction

Material	stainless steel fine turned
Is weight set	yes
Number of weights in set	11
Type of case	in wooden case

Technical data - Measuring system

Nominal value	1 g - 200 g
OIML class	F2
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	160×98×58 mm
Gross weight	1.2 kg
Net weight	1.2 kg

Shipping method Parcel service

Packaging component -
by weight - plastics 0 g

Delivery time 1 d

Optional Services

DAkkS Calibration 962-417

FUNCTIONS

Standard



Option

DAkkS

SERVICES

Model	Description
962-417	DAkkS Calibration KERN 962-417
962-417R	DAkkS Recalibration KERN 962-417R
921-417	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-417
920-712	Marking (Last digits of serial number) KERN 920-712

SPARE PARTS

Model	Description
354-440-200	Box of Weights KERN 354-440-200
315-242	Tweezers KERN 315-242