

KERN ORA 32BA

Refractive index measurement for laboratories and the industry



Article codes, product groups & Status information

Model article-code	ORA 32BA
Model-Serie	ORA-B
Assortment	Microscopes & Optical Metrology
Product Type	Refractometers
Product group	Analogue refractometers
Product usage type	Main article
EAN code (Model)	4045761487130
Predecessor Model	ORA 32BA-2.024a
Planned start date	2025-1HJ
Customs tariff number	90275000
Country of Origin	CN

Technical data - Primary

Scales	Brix
Scale Range	0 % - 32 %
ATC	Yes

Technical data - Construction

Type	Handheld refractometer
Display scale	Analogue
Dimensions (W×D×H)	178×40×40 mm

Technical data - Measuring system

Measuring method	Light refraction
Measuring Temp min	10 °C
Measuring Temp max	30 °C

Functions

Standard

Readout 0,2 %

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	215×94×65 mm
Gross weight	0,45 kg
Net weight	0,20 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	109 g
Packaging component - by weight - plastics	109 g
Delivery time	1 d

Optional Services

ISO Calibration	961-290
-----------------	---------

KERN ORA 32BA

Refractive index measurement for laboratories and the industry



Option

