

KERN 357-13

Individual weight, OIML Class M2, 5 kg, Knob, stainless steel fine turned



Article codes, product groups & Status information

Model article-code	357-13
Model-Serie	357-13
Assortment	KERN
Product Type	Test weights
Product group	OIML M2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761289048
Planned start date	2022-1HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	5 kg
Tolerance (OIML)	800 mg
Case aluminium	317-130-600
Case plastic	347-130-400
Case wood	337-130-200

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	5 kg
OIML class	M2
OIML conform	Yes
Tolerance (OIML)	800 mg

Functions

Standard

Density (kg/m ³)	7.950
Volume (cm ³)	595,2381

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	100×100×165 mm
Gross weight	5,0 kg
Net weight	5,0 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-643
-------------------	---------

KERN 357-13

Individual weight, OIML Class M2, 5 kg, Knob, stainless steel fine turned



Option

