

PRODUCT BROCHURE

OPTICS OZB-A5211

Date: 30.04.2025

ENGLISH



Individuality, variety and flexible working through our modular construction system ? PREMIUM-Universal stands

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

Stereomicroscope Stand OZB-A5211

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OZB-A5211
Model-Serie	OZB-UP
Assortment	Optics
Product Type	Microscopes
Product group	Stereo microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761179271
Customs tariff number	90119000

Gross weight	7 kg
Net weight	6 kg
Shipping method	Parcel service
Packaging component - by weight - paper	0 g
Packaging component - by weight - cardboard	699.75 g
Packaging component - by weight - plastics	23.056 g
Packaging component - by weight - styrofoam	429.994 g
Packaging component - by weight - foam	0 g
Delivery time	1 d

ENGLISH

Technical data - Construction

Stand type	universal (telescopic arm with clamp)
Holder diameter	76 mm
Dimensions (W×D×H)	80×614×515 mm

FUNCTIONS

Standard



Technical data - Environmental conditions

Storage Temp Min	-5 °C
Stotage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	680×330×260 mm
------------------------------	----------------

ACCESSORIES

Model	Description
OZB-A5301	Stereo Microscope Modular System – holders KERN OZB-A5301
OZB-A5306	Stereo Microscope Modular System – holders KERN OZB-A5306

SPARE PARTS

Model	Description
OZB-A5806	Replacement lever OZB-A5806