

KERN 315-207-204

Safety Set 10 g (F1); 100 g (F1)



Article codes, product groups & Status information

Model article-code	315-207-204
Model-Serie	315
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761245471
Planned start date	2019-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	10 g & 100 g
Case aluminium	Zet315-000-000
Weight set case	In aluminium case

Technical data - Construction

Material	Stainless steel
Is weight set	Yes
Number of weights in set	2
Type of case	In aluminium case

Technical data - Measuring system

Nominal value	10 g & 100 g
OIML class	F1
OIML conform	Yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	150×115×70 mm
Gross weight	1.00 kg
Net weight	0.50 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	500 g
Delivery time	1 d

Technical data - Verification

Verification possible	Yes
-----------------------	-----

Optional Services

DAkKS Calibration	962-400-400
-------------------	-------------

Functions

KERN 315-207-204

Safety Set 10 g (F1); 100 g (F1)



Standard

