

KERN YKT-01

High-quality printer with large LCD display and numerous functions



Article codes, product groups & Status information

Model article-code	YKT-01
Model-Serie	YKT
Assortment	KERN
Product Type	General weighing accessories
Product group	Accessories (printer, etc.)
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761484849
Predecessor Model	911-013; YKN-01; YKS-01; YKT-01-2023a
Planned start date	2025-2HJ
Customs tariff number	84239010
CE Label	Yes

Technical data - Construction

Overall dimensions mounted (W×D×H)	114.6×188.2×86.6 mm
------------------------------------	---------------------

Technical data - Power supply

Mains adapter/adapter included	AUS;CH;EURO;UK;US
Input Voltage	100 V - 240 V AC, 50 / 60 Hz
Mains adapter	Power supply unit
Interfaces	RS-232 standard

Technical data - Environmental conditions

Maximum humidity	80 %
------------------	------

Maximum operating temperature	45 °C
-------------------------------	-------

Minimum ambient temperature	5 °C
-----------------------------	------

Minimum humidity	10 %
------------------	------

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	265×200×180 mm
------------------------------	----------------

Gross weight	1.8 kg
--------------	--------

Net weight	0.65 kg
------------	---------

Shipping method	Parcel service
-----------------	----------------

Packaging component - by weight - cardboard	521.5 g
---	---------

Packaging component - by weight - styrofoam	521.5 g
---	---------

Delivery time	1 d
---------------	-----

Functions

Standard

KERN YKT-01

High-quality printer with large LCD display and numerous functions

