

PRODUCT BROCHURE

KERN 347-036

Date: 10.03.2026

ENGLISH



Individual weight, OIML Class M1, 5 g, Hook weight, stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- Individual weight, OIML Class M1, 5 g, Hook weight, stainless steelfine turned (OIML)

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	347-036
Model-Serie	347-036
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761256170
Predecessor Model	347-436
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	5 g
Tolerance (OIML)	1,6 mg
Case plastic	347-030-400

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	5 g
OIML class	M1
OIML conform	yes
Tolerance (OIML)	1,6 mg
Density (kg/m ³)	7950
Volume (cm ³)	0,5952

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	8×8×22 mm
Gross weight	0.05 kg
Net weight	0.05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-633
-------------------	---------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>347-030-400</u>	Plastic Weight Box KERN 347-030-400
<u>315-242</u>	Tweezers KERN 315-242

SERVICES

Model	Description
<u>962-633</u>	DAkkS Calibration KERN 962-633
<u>962-633R</u>	DAkkS Recalibration KERN 962-633R