

PRODUCT BROCHURE

KERN 347-147-100

Date: 10.03.2026

ENGLISH



Supporting bar, OIML Class M1, 10 kg, stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- for fixing slotted weights
- A precise curvature at the lower end of the beam bar ensures that the slotted weights rest securely. This design prevents the weights from slipping and ensures stability – even in slight inclined positions
- Own weight beam bar 10 kg

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	347-147-100
Model-Serie	347-147-100
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 1366855
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	10 kg
Tolerance (OIML)	500 mg

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	10 kg
OIML class	M1
OIML conform	yes
Tolerance (OIML)	500 mg
Density (kg/m ³)	8000
Volume (cm ³)	1190,4762

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	178×178×901.9 mm
Gross weight	10 kg
Net weight	10 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-644
-------------------	---------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>317-290</u>	Gloves KERN 317-290
<u>347-037</u>	Test Weight KERN 347-037
<u>347-047</u>	Test Weight KERN 347-047
<u>347-057</u>	Test Weight KERN 347-057
<u>347-067</u>	Test Weight KERN 347-067
<u>347-077</u>	Test Weight KERN 347-077
<u>347-087</u>	Test Weight KERN 347-087
<u>347-097</u>	Test Weight KERN 347-097
<u>347-117</u>	Test Weight KERN 347-117
<u>347-127</u>	Test Weight KERN 347-127
<u>347-137</u>	Test Weight KERN 347-137
<u>347-147</u>	Test Weight KERN 347-147

SERVICES

Model	Description
<u>962-644</u>	DAkkS Calibration KERN 962-644
<u>962-644R</u>	DAkkS Recalibration KERN 962-644R