

**KERN
SEQ-01**

ENGLISH

ation value (e)	0.001 g	N/A	N/A
	0.01 g	N/A	N/A
ected accessories:	Printer YKG-01		
ing instructions managed by:	Process owner		
Installation:	29-08-2023		
putting out of service:	N/A		
contract number:	N/A		
case telephone number:	N/A		
ment qualification executed by:	QUALIFICATION OFFICER 		
Günther			
of qualification officer)	(signature)		
ment qualification reviewed by:			
Doe			
of measuring and test equipment officer)	(signature)		

Equipment Qualification made easy – with the KERN Equipment Qualification Software (EQS)

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more information

DESCRIPTION

- Equipment qualification is far more than just an administrative obligation. Equipment qualification for balances checks whether the balance is suitable for the intended process with the technology being used.
- Equipment qualification can be carried out with a paper-based system, but the amount of time and paperwork required to do this is not insignificant. This process is significantly easier and quicker using the KERN EQS software. With software support there is the option of entering master data e.g. about the balance or the user just once and then copying this to all relevant places in the program. Even better: Important master data about the balances such as, for example maximum load or readout, can be taken in full from existing KERN systems by the KERN EQS software.

Benefits

- One-time input of master data: Only input the master data such as balance information or user details once. Our software will distribute this information automatically to all relevant places within a qualification procedure
- Automatic data transfer: Important balance data such as maximum load and readout will be taken straight from your existing KERN systems. No manual entry required, which eliminates possible sources of error
- Efficient qualification: Our software will guide you through the qualification process, by providing required fields and optional fields. This ensures that all important information is available, and that the processing time is kept to a minimum due to one-time input
- Signatures: These are entered and stored within the KERN EQS software and can then be printed out.
- User roles: Clear distinction between test users and review users to make sure there is no non-compliance with the double-check principle
- Deviation list: Deviations or errors immediately go onto a task list, which can be worked through

Licensing

- License period: 2 years – you can purchase a license extension in combination with a KERN qualification training course.
- Scope of the software license: Can be installed on up to 4 computers (PCs, tablets) or just one (trained) person. Can be used on as many balances as required

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	SEQ-01
--------------------	--------

Model-Serie	SEQ
Assortment	KERN
Product Type	Software

Product group	Software (balances), other
Product usage type	Main article
EAN code (Model)	4045761361720
Customs tariff number	85234910
Country of Origin	DE

Technical data - Packing & Shipping

Packaging component - by weight - paper	100 g
Packaging component - by weight - cardboard	0 g
Packaging component - by weight - styrofoam	0 g
Packaging component - by weight - glass	0 g
Packaging component - by weight - wood	0 g
Packaging component - by weight - foam	0 g
Packaging component - by weight - metal	0 g
Delivery time	1 d