

KERN 357-01

Individual weight, OIML Class M2, 1 g, Knob, stainless steel fine turned



Article codes, product groups & Status information

Model article-code	357-01
Model-Serie	357-01
Assortment	KERN
Product Type	Test weights
Product group	OIML M2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761288836
Planned start date	2022-1HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 g
Tolerance (OIML)	3 mg
Case aluminium	317-010-600
Case plastic	347-030-400
Case wood	337-010-200

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	1 g
OIML class	M2
OIML conform	Yes
Tolerance (OIML)	3 mg

Functions

Standard

Density (kg/m³) 7.950

Volume (cm³) 0,119

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	24×24×51 mm
Gross weight	0,05 kg
Net weight	0,05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkKS Calibration 962-631

KERN 357-01

Individual weight, OIML Class M2, 1 g, Knob, stainless steel fine turned



Option

