

PRODUCT BROCHURE

KERN 343-036

Date: 10.03.2026

ENGLISH



Weight set, M1, 1 mg - 100 g, Knob, stainless steel fine turned, in aluminium case

CONTACT


KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



DESCRIPTION

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	343-036
Model-Serie	343
Assortment	KERN
Product Type	Test weights
Product group	OIML M 1 Test weights
Product usage type	Main article
EAN code (Model)	4045761224636
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 mg - 100 g
---------------	--------------

Case aluminium	313-030-600
Weight set case	in aluminium case

Technical data - Construction

Material	stainless steel fine turned
Is weight set	yes
Type of case	in aluminium case

Technical data - Measuring system

Nominal value	1 mg - 100 g
OIML class	M1
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	165×125×75 mm
Gross weight	0.90 kg
Net weight	0.70 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	240 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-602
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
<u>962-602</u>	DAkkS Calibration KERN 962-602
<u>962-602R</u>	DAkkS Recalibration KERN 962-602R
<u>920-703</u>	Marking (Last digits of serial number) KERN 920-703

SPARE PARTS

Model	Description
<u>313-030-600</u>	Aluminium Protective Box KERN 313-030-600