


KERN 338-01

Milligram weights, F2, 1 mg, Platelet, stainless steel polished



Article codes, product groups & Status information

Model article-code	338-01
Model-Serie	338
Assortment	KERN
Product Type	Test weights
Product group	OIML F2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761240223
Predecessor Model	338-1-219a
Planned start date	2020-1HJ
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary

Nominal value	1 mg
Tolerance (OIML)	0,06 mg
Case aluminium	317-009-600
Case plastic	347-009-400
Case wood	338-090-200

Technical data - Construction

Dimensions (W×D×H)	3,5×3,5×3 mm
Material	Stainless steel polished

Technical data - Measuring system

Nominal value	1 mg
OIML class	F2

Functions

Standard

OIML conform	Yes
Tolerance (OIML)	0,06 mg
Density (kg/m ³)	7.960
Volume (cm ³)	0,0004

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	3,5×3,5×3 mm
Gross weight	0,05 kg
Net weight	0,05 kg
Shipping method	Parcel service
Delivery time	1 d

Optional Services

DAkkS Calibration	962-451
-------------------	---------

KERN 338-01

Milligram weights, F2, 1 mg, Platelet, stainless steel polished



Option

