

PRODUCT BROCHURE

KERN 367-06

Date: 10.03.2026

ENGLISH



Individual weight, OIML Class M3, 50 g, Knob, stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	367-06
Model-Serie	367-06
Assortment	KERN
Product Type	Test weights
Product group	OIML M3 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761257313
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	50 g
Tolerance (OIML)	30 mg
Case aluminium	317-060-600
Case plastic	347-070-400
Case wood	337-060-200

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	50 g
OIML class	M3

OIML conform	yes
Tolerance (OIML)	30 mg
Density (kg/m ³)	7950
Volume (cm ³)	5,9524

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	18×18×28 mm
Gross weight	0.05 kg
Net weight	0.05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	no
-----------------------	----

Optional Services

DAkkS Calibration	962-636
-------------------	---------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>337-060-200</u>	Wooden Box for Set of Weights KERN 337-060-200
<u>347-070-400</u>	Plastic Weight Box KERN 347-070-400
<u>317-060-600</u>	Aluminium Protective Box KERN 317-060-600
<u>315-242</u>	Tweezers KERN 315-242

SERVICES

Model	Description
<u>962-636</u>	DAkkS Calibration KERN 962-636
<u>962-636R</u>	DAkkS Recalibration KERN 962-636R