


# KERN 347-022-005

Weight container for test weights, calibrated to OIML class M1



## Article codes, product groups & Status information

|                       |                      |
|-----------------------|----------------------|
| Model article-code    | 347-022-005          |
| Model-Serie           | 347                  |
| Assortment            | KERN                 |
| Product Type          | Test weights         |
| Product group         | OIML M1 Test weights |
| Product usage type    | Main article         |
| EAN code (Model)      | 4045761225015        |
| Planned start date    | 2019-1HJ             |
| Customs tariff number | 84239090             |

|                                            |                |
|--------------------------------------------|----------------|
| Net weight                                 | 20 kg          |
| Shipping method                            | Parcel service |
| Packaging component - by weight - plastics | 0 g            |
| Delivery time                              | 1 d            |

## Optional Services

|                   |         |
|-------------------|---------|
| DAkkS Calibration | 962-645 |
|-------------------|---------|

## Technical data - Primary

|                  |       |
|------------------|-------|
| Nominal value    | 20 kg |
| Tolerance (OIML) | 1 g   |

## Technical data - Construction

|                    |                 |
|--------------------|-----------------|
| Dimensions (W×D×H) | 1000×500×500 mm |
| Material           | Stainless steel |

## Technical data - Measuring system

|                  |       |
|------------------|-------|
| Nominal value    | 20 kg |
| OIML class       | M1    |
| OIML conform     | Yes   |
| Tolerance (OIML) | 1 mg  |

## Technical data - Packing & Shipping

|                              |                 |
|------------------------------|-----------------|
| Dimensions packaging (W×D×H) | 1000×500×500 mm |
| Gross weight                 | 20 kg           |

## Functions

# KERN 347-022-005

Weight container for test weights, calibrated to OIML class M1



Standard

