

KERN 326-092

Weight set, F1, 1 g - 10 kg, Eco Shape, stainless steel polished, in wooden case



Article codes, product groups & Status information

Model article-code	326-092
Model-Serie	326
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main article
EAN code (Model)	4045761224315
Planned start date	2019-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	1 g - 10 kg
Case wood	313-090-100
Weight set case	In wooden case

Technical data - Construction

Material	Stainless steel polished
Is weight set	Yes
Type of case	In wooden case

Technical data - Measuring system

Nominal value	1 g - 10 kg
OIML class	F1
OIML conform	Yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	540×200×240 mm
------------------------------	----------------

Gross weight	34 kg
Net weight	32 kg
Shipping method	Freight
Packaging component - by weight - plastics	1890 g
Delivery time	2 d

Optional Services

DAkkS Calibration	962-422
-------------------	---------

Functions

KERN 326-092

Weight set, F1, 1 g - 10 kg, Eco Shape, stainless steel polished, in wooden case



Standard

