

KERN OZB-A4516

Double gooseneck illumination for best illumination of your sample



Article codes, product groups & Status information

Model article-code	OZB-A4516
Model-Serie	OZB-IF
Assortment	Optics
Product Type	Microscopes
Product group	External Illumination
Product usage type	Main article
EAN code (Model)	4045761217720
Predecessor Model	OZB-A4512
Planned start date	2019-1HJ
Customs tariff number	94054239

Technical data - Primary

Illumination type	LED
Illumination intensity	20W

Technical data - Construction

Color tempertaure	6400 K
Dimmable	Yes
Fibre length	540 mm
Dimensions	100×200×150 mm

Technical data - Power supply

Protection Class	Class III - Low voltage
------------------	-------------------------

Technical data - Environmental conditions

Storage Temp Min	-5 °C
------------------	-------

Functions

Standard

Stotage Temp Max	40 °C
------------------	-------

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	300×250×250 mm
Gross weight	4 kg
Net weight	3.8 kg
Shipping method	Parcel service
Packaging component - by weight - paper	0 g
Packaging component - by weight - cardboard	207.576 g
Packaging component - by weight - plastics	1.41932 g
Packaging component - by weight - styrofoam	145.835 g
Packaging component - by weight - foam	0 g
Delivery time	1 d

Technical data - Verification

WEEE	Yes
------	-----

KERN OZB-A4516

Double gooseneck illumination for best illumination of your sample

