

KERN 357-09

Individual weight, OIML Class M2, 500 g, Knob, stainless steel fine turned



Article codes, product groups & Status information

Model article-code	357-09
Model-Serie	357-09
Assortment	KERN
Product Type	Test weights
Product group	OIML M2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761289017
Planned start date	2022-1HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	500 g
Tolerance (OIML)	80 mg
Case aluminium	317-090-600
Case plastic	347-090-400
Case wood	337-090-200

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	500 g
OIML class	M2
OIML conform	Yes
Tolerance (OIML)	80 mg

Functions

Standard

Density (kg/m ³)	7.950
Volume (cm ³)	59,5238

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	57×57×74 mm
Gross weight	0,50 kg
Net weight	0,50 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkKS Calibration	962-639
-------------------	---------

KERN 357-09

Individual weight, OIML Class M2, 500 g, Knob, stainless steel fine turned



Option

