

PRODUCT BROCHURE

KERN
328-39

Date: 02.05.2026

ENGLISH



Milligram weights, F1, 500 mg, Wire, stainless steel polished

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	328-39
Model-Serie	328
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 1478664
Predecessor Model	328-39-2024a
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary

Nominal value	500 mg
Tolerance (OIML)	0,08 mg
Case aluminium	317-009-600
Case plastic	347-009-400
Case wood	338-090-200

Technical data - Construction

Dimensions (W×D×H)	15×15×4 mm
Material	stainless steel polished

Technical data - Measuring system

Nominal value	500 mg
---------------	--------

OIML class F1

OIML conform yes

Tolerance (OIML) 0,08 mg

Density (kg/m³) 8000

Volume (cm³) 0,0625

Technical data - Packing & Shipping

Dimensions packaging (W×D×H) 15×15×4 mm

Gross weight 0.05 kg

Net weight 0.05 kg

Shipping method Parcel service

Packaging component - by weight - plastics 0 g

Delivery time 1 d

Optional Services

DAkkS Calibration 962-459

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>338-090-200</u>	Wooden Box for Set of Weights KERN 338-090-200
<u>347-009-400</u>	Plastic Weight Box KERN 347-009-400
<u>317-009-600</u>	Aluminium Protective Box KERN 317-009-600
<u>315-242</u>	Tweezers KERN 315-242
<u>318-270</u>	Dusting Brush KERN 318-270

SERVICES

Model	Description
<u>962-459</u>	DAkkS Calibration KERN 962-459
<u>962-459R</u>	DAkkS Recalibration KERN 962-459R

RELATED ARTICLES DATA

Model name	Nominal value	Tolerance (OIML)	Case aluminium	Case plastic	Case wood	Weight set case
<u>328-226</u>	1 mg - 500 mg	-	313-010-600	-	-	in aluminium case
<u>328-01</u>	1 mg	0,02 mg	317-009-600	347-009-400	338-090-200	-
<u>328-02</u>	2 mg	0,02 mg	317-009-600	347-009-400	338-090-200	-
<u>328-03</u>	5 mg	0,02 mg	317-009-600	347-009-400	338-090-200	-
<u>328-04</u>	10 mg	0,025 mg	317-009-600	347-009-400	338-090-200	-
<u>328-05</u>	20 mg	0,03 mg	317-009-600	347-009-400	338-090-200	-
<u>328-06</u>	50 mg	0,04 mg	317-009-600	347-009-400	338-090-200	-
<u>328-07</u>	100 mg	0,05 mg	317-009-600	347-009-400	338-090-200	-
<u>328-08</u>	200 mg	0,06 mg	317-009-600	347-009-400	338-090-200	-
<u>328-09</u>	500 mg	0,08 mg	317-009-600	347-009-400	338-090-200	-
<u>328-22</u>	1 mg - 500 mg	-	-	328-222-420	-	in plastic case
<u>328-39</u>	500 mg	0,08 mg	317-009-600	347-009-400	338-090-200	-
<u>328-38</u>	200 mg	0,06 mg	317-009-600	347-009-400	338-090-200	-
<u>328-37</u>	100 mg	0,05 mg	317-009-600	347-009-400	338-090-200	-
<u>328-36</u>	50 mg	0,04 mg	317-009-600	347-009-400	338-090-200	-
<u>328-35</u>	20 mg	0,03 mg	317-009-600	347-009-400	338-090-200	-

Model name	Nominal value	Tolerance (OIML)	Case aluminium	Case plastic	Case wood	Weight set case
<u>328-34</u>	10 mg	0,025 mg	317-009-600	347-009-400	338-090-200	-
<u>328-33</u>	5 mg	0,02 mg	317-009-600	347-009-400	338-090-200	-
<u>328-32</u>	2 mg	0,02 mg	317-009-600	347-009-400	338-090-200	-
<u>328-31</u>	1 mg	0,02 mg	317-009-600	347-009-400	338-090-200	-