

PRODUCT BROCHURE

KERN 304-084

Date: 10.03.2026

ENGLISH



Weight set, E1, 1 g - 5 kg, Knob, stainless steel polished, in plastic case

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	304-084
Model-Serie	304
Assortment	KERN
Product Type	Test weights
Product group	OIML E1 Test weights
Product usage type	Main article
EAN code (Model)	4045761291744
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 g - 5 kg
Case plastic	313-082-400
Weight set case	in plastic case

Technical data - Construction

Material	stainless steel polished
Is weight set	yes
Number of weights in set	16
Type of case	in plastic case

Technical data - Measuring system

Nominal value	1 g - 5 kg
OIML class	E1
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	350×320×180 mm
Gross weight	13 kg
Net weight	13 kg
Shipping method	Parcel service

Packaging component - 0 g
by weight - plastics

Delivery time 1 d

Optional Services

DAkkS Calibration 963-221

FUNCTIONS

Standard



Option



SERVICES

Model	Description
963-221	DAkkS Calibration KERN 963-221
962-221R	DAkkS Recalibration KERN 962-221R
920-716	Marking (Last digits of serial number) KERN 920-716

SPARE PARTS

Model	Description
313-082-400	Plastic Weight Box KERN 313-082-400
315-242	Tweezers KERN 315-242