


KERN 328-03

Milligram weights, F1, 5 mg, Platelet, stainless steel polished



Article codes, product groups & Status information

Model article-code	328-03
Model-Serie	328
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761240339
Predecessor Model	328-3-2.019a
Planned start date	2020-1HJ
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary

Nominal value	5 mg
Tolerance (OIML)	0,02 mg
Case aluminium	317-009-600
Case plastic	347-009-400
Case wood	338-090-200

Technical data - Construction

Dimensions (W×D×H)	6×6×3 mm
Material	Stainless steel polished

Technical data - Measuring system

Nominal value	5 mg
OIML class	F1

Functions

Standard

OIML conform	Yes
Tolerance (OIML)	0,02 mg
Density (kg/m ³)	7.960
Volume (cm ³)	0,0019

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	6×6×3 mm
Gross weight	0,05 kg
Net weight	0,05 kg
Shipping method	Parcel service
Delivery time	1 d

Optional Services

DAkKS Calibration	962-453
-------------------	---------

KERN 328-03

Milligram weights, F1, 5 mg, Platelet, stainless steel polished



Option

