

## PRODUCT BROCHURE

---

# OPTICS OZB-A5223

Date: 22.09.2025

ENGLISH



Individuality, variety and flexible working through our modular construction system ? PREMIUM-Universal stands

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## PRODUCT SHOWCASE



## DESCRIPTION

Stereomicroscope Stand OZB-A5223

## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                                 |
|-----------------------|---------------------------------|
| Model article-code    | OZB-A5223                       |
| Model-Serie           | OZB-UP                          |
| Assortment            | Optics                          |
| Product Type          | Microscopes                     |
| Product group         | Stereo microscopes              |
| Product usage type    | Spare part;Accessories/Services |
| EAN code (Model)      | 404576 12 17935                 |
| Customs tariff number | 90 11 9000                      |
| CE Label              | yes                             |

### Technical data - Construction

|            |                                                            |
|------------|------------------------------------------------------------|
| Stand type | universal (ball-bearing double telescopic arm with screws) |
|------------|------------------------------------------------------------|

|                    |                |
|--------------------|----------------|
| Stand area         | 150×150 mm     |
| Holder diameter    | 76 mm          |
| Dimensions (W×D×H) | 150×545×515 mm |

### Technical data - Environmental conditions

|                  |       |
|------------------|-------|
| Storage Temp Min | -5 °C |
| Storage Temp Max | 40 °C |

### Technical data - Packing & Shipping

|                                         |                |
|-----------------------------------------|----------------|
| Dimensions packaging (W×D×H)            | 680×330×260 mm |
| Gross weight                            | 11 kg          |
| Net weight                              | 9 kg           |
| Shipping method                         | Parcel service |
| Packaging component - by weight - paper | 0 g            |

Packaging component - 1140.53 g  
by weight - cardboard

Packaging component - 22.9164 g  
by weight - plastics

Packaging component - 599.352 g  
by weight - styrofoam

Packaging component - 0 g  
by weight - foam

Delivery time 1 d

## FUNCTIONS

Standard



## ACCESSORIES

| Model     | Description                                               |
|-----------|-----------------------------------------------------------|
| OZB-A5301 | Stereo Microscope Modular System – holders KERN OZB-A5301 |
| OZB-A5306 | Stereo Microscope Modular System – holders KERN OZB-A5306 |

## SPARE PARTS

| Model     | Description                                |
|-----------|--------------------------------------------|
| OZB-A5806 | Replacement lever OZB-A5806                |
| OZB-A5292 | Clamping plates with screws KERN OZB-A5292 |