


PRODUCT BROCHURE

KERN YKA-48

Date: 16.03.2025

ENGLISH



Power adapter (CH, EURO, UK)

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

Power adapter YKA-48 (5 V 1000 mA | 100 V - 240 V, 50 / 60 Hz | EURO, UK, CH)

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	YKA-48
Model-Serie	YKA
Assortment	KERN
Product Type	General weighing accessories
Product group	Accessories (printer, etc.)
Product usage type	Spare part;Accessories/Services
EAN code (Model)	404576 1357 18 1
Customs tariff number	85044095
CE Label	yes

Packaging component -
by weight - cardboard 20 g

Packaging component -
by weight - plastics 0 g

Delivery time 1 d

FUNCTIONS

Standard



ENGLISH

Technical data - Power supply

Mains adapter/adapter included	CH;EURO;UK
Input Voltage	100 V - 240 V, 50 / 60 Hz

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	90×80×50 mm
Gross weight	0.15 kg
Net weight	0.10 kg
Shipping method	Parcel service