

PRODUCT BROCHURE

KERN YKA-10

Date: 16.03.2025

ENGLISH



Mains adapter 100 V - 240 V, 50 / 60 Hz, standard EURO/US/AUS/UK/CH

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

Power adapter YKA-10 (15 V, 2.1 A | 100 V - 240 V, 50 / 60 Hz | EURO, UK, US)

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	YKA-10
Model-Serie	YKA
Assortment	KERN
Product Type	General weighing accessories
Product group	Accessories (printer, etc.)
Product usage type	Spare part;Accessories/Services
EAN code (Model)	404576 1174344
Predecessor Model	PLS-A13
Customs tariff number	85044095
CE Label	yes

Packaging component - by weight - cardboard	72.5 g
---	--------

Packaging component - by weight - plastics	72.5 g
--	--------

Packaging component - by weight - glass	0 g
---	-----

Packaging component - by weight - wood	0 g
--	-----

Delivery time	1 d
---------------	-----

FUNCTIONS

Standard



ENGLISH

Technical data - Power supply

Mains adapter/adapter included	EURO;UK;US
Input Voltage	100 V - 240 V, 50 / 60 Hz

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	70×70×70 mm
Gross weight	0.90 kg
Net weight	0.75 kg
Shipping method	Parcel service