

PRODUCT BROCHURE

KERN 354-07

Date: 10.03.2026

ENGLISH



Weight set, M2, 1 g - 2 kg, Knob, stainless steel fine turned, in wooden case

CONTACT


KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	354-07
Model-Serie	354
Assortment	KERN
Product Type	Test weights
Product group	OIML M2 Test weights
Product usage type	Main article
EAN code (Model)	404576 1292390
Predecessor Model	354-47
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 g - 2 kg
Case wood	353-470-200
Weight set case	in wooden case

Technical data - Construction

Material	stainless steel fine turned
Is weight set	yes
Number of weights in set	15
Type of case	in wooden case

Technical data - Measuring system

Nominal value	1 g - 2 kg
OIML class	M2
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	305×180×140 mm
Gross weight	10 kg
Net weight	10 kg

Shipping method Parcel service

Packaging component -
by weight - plastics 0 g

Delivery time 1 d

Optional Services

DAkkS Calibration 962-620

FUNCTIONS

Standard



Option

DAkkS

SERVICES

Model	Description
962-620	DAkkS Calibration KERN 962-620
962-620R	DAkkS Recalibration KERN 962-620R
920-715	Marking (Last digits of serial number) KERN 920-715

SPARE PARTS

Model	Description
353-470-200	Box of Weights KERN 353-470-200
315-242	Tweezers KERN 315-242