

KERN 316-11

Individual weight, OIML Class E2, 1 kg, Compact form, stainless steel polished



Article codes, product groups & Status information

Model article-code	316-11
Model-Serie	316
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761094925
Planned start date	2004-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	1 kg
Tolerance (OIML)	1,6 mg
Case aluminium	317-110-600
Case plastic	317-110-400
Case wood	317-110-100

Technical data - Construction

Material	Stainless steel polished
----------	--------------------------

Technical data - Measuring system

Nominal value	1 kg
OIML class	E2
OIML conform	Yes
Tolerance (OIML)	1,6 mg
Density (kg/m ³)	7950

Functions

Standard

Volume (cm³) 125

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	48×48×73 mm
Gross weight	1.00 kg
Net weight	1.00 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	No
-----------------------	----

Optional Services

DAkkS Calibration	962-341
-------------------	---------

KERN 316-11

Individual weight, OIML Class E2, 1 kg, Compact form, stainless steel polished

