

PRODUCT BROCHURE

KERN YKI-10

Date: 19.09.2025

ENGLISH



RS-232/Ethernet-Converter

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



DESCRIPTION

- Converter for connecting scales, force measuring instruments, etc. with an RS-232 data interface to an IP-based Ethernet network
- Allows direct integration of measuring devices into most company networks, with automatic configuration using DHCP or statistical IP allocation. The current measurement data can be accessed from the entire network at any time
- The converter supports a wide range of configuration options via the web interface
- Compatible with KERN BalanceConnection software
- Suitable for wall mounting or free-standing
- Included with delivery: RS-232/Ethernet converter with 1 × RS-232 port, mains adapter (EU)

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	YKI-10
Model-Serie	YKI
Assortment	KERN
Product Type	General weighing accessories
Product group	Accessories (printer, etc.)
Product usage type	Accessories/DL

EAN code (Model)	4045761483361
------------------	---------------

Predecessor Model	YKI-01
-------------------	--------

Customs tariff number	84239010
-----------------------	----------

CE Label	yes
----------	-----

Technical data - Power supply


Mains adapter	Power supply unit
---------------	-------------------

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	230×165×50 mm
Gross weight	0.35 kg
Net weight	0.25 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	130 g
Delivery time	1 d

FUNCTIONS

Standard



SPARE PARTS

Model	Description
YKA-11	Power supply KERN YKA-11