

PRODUCT BROCHURE

OPTICS OZB-A1201

Date: 30.04.2025

ENGLISH



Individuality, variety and flexible working through our modular construction system ? ECO-Universal stands

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

Stereomicroscope Stand OZB-A1201

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OZB-A1201
Model-Serie	OZB-UE
Assortment	Optics
Product Type	Microscopes
Product group	Stereo microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761218284
Customs tariff number	90119000
CE Label	yes

Technical data - Construction

Stand type	universal (telescopic arm with plate)
Stand area	230×230 mm
Holder diameter	76 mm
Dimensions (W×D×H)	410×230×420 mm

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Stotage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	490×420×290 mm
------------------------------	----------------

Gross weight	16 kg
--------------	-------

Net weight	12 kg
------------	-------

Shipping method	Parcel service
-----------------	----------------

Packaging component - by weight - cardboard	2 kg
---	------

Packaging component - by weight - plastics	2 kg
--	------

Delivery time	1 d
---------------	-----

FUNCTIONS

Standard



ENGLISH

SPARE PARTS

Model	Description
OZB-A1904	Replacement screw OZB-A1904
OZB-A1905	Holder connecting rod OZB-A1905