

KERN 367-08

Individual weight, OIML Class M3, 200 g, Knob, stainless steel fine turned



Article codes, product groups & Status information

Model article-code	367-08
Model-Serie	367-08
Assortment	KERN
Product Type	Test weights
Product group	OIML M3 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761257276
Planned start date	2022-1HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	200 g
Tolerance (OIML)	100 mg
Case aluminium	317-080-600
Case plastic	347-080-400
Case wood	337-080-200

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	200 g
OIML class	M3
OIML conform	Yes
Tolerance (OIML)	100 mg

Functions

Standard

Density (kg/m ³)	7.950
Volume (cm ³)	23,8095

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	28×28×47 mm
Gross weight	0,20 kg
Net weight	0,20 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	No
-----------------------	----

Optional Services

DAkkS Calibration	962-638
-------------------	---------

KERN 367-08

Individual weight, OIML Class M3, 200 g, Knob, stainless steel fine turned



Option

