


PRODUCT BROCHURE

OPTICS OZB-A5222

Date: 30.04.2025

ENGLISH



Individuality, variety and flexible working through our modular construction system ? PREMIUM-Universal stands

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



DESCRIPTION

Stereomicroscope Stand OZB-A5222

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OZB-A5222
Model-Serie	OZB-UP
Assortment	Optics
Product Type	Microscopes
Product group	Stereo microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	404576 12 17928
Customs tariff number	90 11 9000
CE Label	yes

Technical data - Construction

Stand type	universal (jointed arm with screws)
------------	-------------------------------------

Stand area	150×150 mm
Holder diameter	76 mm
Dimensions (W×D×H)	150×553×515 mm

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	680×330×260 mm
Gross weight	10 kg
Net weight	9 kg
Shipping method	Parcel service
Packaging component - by weight - paper	0 g

Packaging component - 1 140.28 g
by weight - cardboard

Packaging component - 23.2372 g
by weight - plastics

Packaging component - 496.28 g
by weight - styrofoam

Packaging component - 0 g
by weight - foam

Delivery time 1 d

FUNCTIONS

Standard



ACCESSORIES

Model	Description
OZB-A5301	Stereo Microscope Modular System – holders KERN OZB-A5301
OZB-A5306	Stereo Microscope Modular System – holders KERN OZB-A5306

SPARE PARTS

Model	Description
OZB-A5806	Replacement lever OZB-A5806
OZB-A5292	Clamping plates with screws KERN OZB-A5292