

OPTICS ORL-A2002

Date: 28.01.2025

ENGLISH



USB charger

CONTACT

 Address

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

 Phone and fax

Phone : +49 - [0] 7433 - 9933 - 0
Fax : +49 - [0] 7433 - 9933 - 149

 Internet

Email : info@kern-sohn.com
Web : www.kern-sohn.com

PRODUCT SHOWCASE



DESCRIPTION

USB charger ORL-A2002

SPECIFICATION

Article codes, product groups & Status information

Model article-code	ORL-A2002
Model-Serie	ORL-A
Assortment	Optics
Product Type	Refractometers
Product group	Digital refractometers
Product usage type	Spare part
EAN code (Model)	404576 126345 1
Planned start date	2022-1HJ
Customs tariff number	90279000

Technical data - Construction

Dimensions	45×45×40 mm
------------	-------------

Technical data - Power supply

Power supply kind	USB including power supply
Input Voltage	100 V - 240V AC 50/60Hz 0.3A
Power Plug	USB socket
Protection Class	Protection class II - protective insulation

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Stotage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	80×70×40 mm
Gross weight	0.10 kg
Net weight	0.05 kg
Shipping method	Parcel service

Packaging component - by weight - cardboard	71.7 g
--	--------

Delivery time	1 d
---------------	-----

Technical data - Verification

CE	yes
----	-----

WEEE	yes
------	-----

RELATED ARTICLES DATA

Model name	Scales	Scale Range	ATC
ORD-A2102			
ORD-A2106			
ORA-A2262			
ORM-A2001			
ORL-A2003			
ORL-A2004			
ORL-A2005			
ORL-A2006			
ORL-A2002			
ORM-A2002			
ORM-A2003			
ORM-A2004			
ORM-A2005			
ORL-A2001			

ENGLISH