


PRODUCT BROCHURE

KERN 338-06

Date: 10.03.2026

ENGLISH



Milligram weights, F2, 50 mg, Platelet, stainless steel polished

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- Milligram weights, F2, 50 mg, Platelet, stainless steel polished (OIML)
- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section
- For the verification you need a box. Please add a box for single weights.

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	338-06
Model-Serie	338
Assortment	KERN
Product Type	Test weights
Product group	OIML F2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 1240278
Predecessor Model	338-06-2019a
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary

Nominal value	50 mg
Tolerance (OIML)	0,12 mg
Case aluminium	317-009-600
Case plastic	347-009-400
Case wood	338-090-200

Technical data - Construction

Dimensions (W×D×H)	7,5×6×3 mm
Material	stainless steel polished

Technical data - Measuring system

Nominal value	50 mg
OIML class	F2
OIML conform	yes
Tolerance (OIML)	0,12 mg
Density (kg/m ³)	7960
Volume (cm ³)	0,0058

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	7.5×6×3 mm
Gross weight	0.05 kg
Net weight	0.05 kg
Shipping method	Parcel service
Delivery time	1 d

Technical data - Verification

Verification possible yes

Optional Services

DAkkS Calibration 962-456

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>338-090-200</u>	Wooden Box for Set of Weights KERN 338-090-200
<u>347-009-400</u>	Plastic Weight Box KERN 347-009-400
<u>317-009-600</u>	Aluminium Protective Box KERN 317-009-600
<u>315-242</u>	Tweezers KERN 315-242
<u>318-270</u>	Dusting Brush KERN 318-270

SERVICES

Model	Description
<u>962-456</u>	DAkkS Calibration KERN 962-456
<u>962-456R</u>	DAkkS Recalibration KERN 962-456R
<u>921-456</u>	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-456