

PRODUCT BROCHURE

KERN 325-034

Date: 10.03.2026

ENGLISH



Weight set, F1, 1 mg - 100 g, Eco Shape, stainless steel polished, in plastic case

CONTACT



KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



ENGLISH

DESCRIPTION

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	325-034
Model-Serie	325
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights

Product usage type	Main article
EAN code (Model)	404576 1224056
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 mg - 100 g
---------------	--------------

Case plastic 313-052-400

Weight set case in plastic case

Option

DAkkS

Technical data - Construction

Material stainless steel polished

Is weight set yes

Type of case in plastic case

Technical data - Measuring system

Nominal value 1 mg - 100 g

OIML class F1

OIML conform yes

Technical data - Packing & Shipping

Dimensions packaging
(W×D×H) 250×210×75 mm

Gross weight 0.60 kg

Net weight 0.55 kg

Shipping method Parcel service

Packaging component -
by weight - plastics 70 g

Delivery time 1 d

Optional Services

DAkkS Calibration 962-402

FUNCTIONS

Standard



SERVICES

Model	Description
<u>962-402</u>	DAkkS Calibration KERN 962-402
<u>962-402R</u>	DAkkS Recalibration KERN 962-402R
<u>921-402</u>	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-402
<u>920-703</u>	Marking (Last digits of serial number) KERN 920-703

SPARE PARTS

Model	Description
<u>313-052-400</u>	Plastic Weight Box KERN 313-052-400