

PRODUCT BROCHURE

KERN 333-036

Date: 10.03.2026

ENGLISH



Weight set, F2, 1 mg - 100 g, Knob, stainless steel fine turned, in aluminium case

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



DESCRIPTION

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	333-036
Model-Serie	333
Assortment	KERN
Product Type	Test weights
Product group	OIML F2 Test weights
Product usage type	Main article
EAN code (Model)	404576 122455 1
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 mg - 100 g
---------------	--------------

Case aluminium 313-030-600

Weight set case in aluminium case

Technical data - Construction

Material	stainless steel fine turned
Is weight set	yes
Type of case	in aluminium case

Technical data - Measuring system

Nominal value	1 mg - 100 g
OIML class	F2
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	165×125×75 mm
Gross weight	0.90 kg
Net weight	0.65 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	250 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-402
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
<u>962-402</u>	DAkkS Calibration KERN 962-402
<u>962-402R</u>	DAkkS Recalibration KERN 962-402R
<u>921-402</u>	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-402
<u>920-703</u>	Marking (Last digits of serial number) KERN 920-703

SPARE PARTS

Model	Description
<u>313-030-600</u>	Aluminium Protective Box KERN 313-030-600