

PRODUCT BROCHURE

KERN YKT-01

Date: 19.09.2025

ENGLISH



High-quality printer with large LCD display and numerous functions

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



ENGLISH

DESCRIPTION

- GLP-compliant printout thanks to durable paper and dot matrix printing process
- Printout of date and time to document the weighings with a time stamp
- Specification of name, batch number and ID for printing user-defined information
- Display and printout of tare, gross and net weights for easy calculation and documentation
- Multiple weighing for storing and processing multiple consecutive weighings
- Recipe function for weighing and documenting various components
- Totalizing function for calculating the total of multiple weighings
- Statistical function for calculating and printing statistical values such as average value or standard deviation
- Printing process: 8 print magnets with dot matrix
- Speed: 1,7 line/s
- Character size (ASCII characters): W×H 1.7×2.6 mm
- Service life: 1 million lines
- Ribbon cartridge: EPSON ERC-22 or ERC-09B
- Operation via mains adapter 230 V/50 Hz (standard EU), as standard

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	YKT-01
--------------------	--------

Model-Serie	YKT
Assortment	KERN
Product Type	General weighing accessories

Product group	Accessories (printer, etc.)
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 1484849
Predecessor Model	9 11-013; YKN-01; YKS-01; YKT-01-2023a
Customs tariff number	842390 10
CE Label	yes

Technical data - Construction

Overall dimensions mounted (W×D×H)	114.6×188.2×86.6 mm
------------------------------------	---------------------

Technical data - Power supply

Mains adapter/adapter included	AUS;CH;EURO;UK;US
Input Voltage	100 V - 240 V AC, 50 / 60 Hz
Mains adapter	Power supply unit

Technical data - Environmental conditions

Maximum humidity	80 %
Maximum operating temperature	45 °C
Minimum ambient temperature	5 °C
Minimum humidity	10 %

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	265×200×180 mm
Gross weight	1.8 kg

Net weight	0.65 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	521.5 g
Packaging component - by weight - styrofoam	521.5 g
Delivery time	1 d

FUNCTIONS

Standard



ACCESSORIES

Model	Description
911-013-010	Paper Rolls KERN 911-013-010
911-013-002	Paper Rolls KERN 911-013-002
YKT-A03	Colour Band KERN YKT-A03
YKT-A05S05	KERN YKT-A05S05 label rolls

SERVICES

Model	Description
970-101	Warranty extension KERN (+2 years) 970-101