

PRODUCT BROCHURE

OPTICS OBB-A1617

Date: 03.04.2025

ENGLISH



Eyepiece: HWF 10 x / Ø 20 mm, with scale 0,1 mm, anti-fungus, high-eye-point, (with scale 0,1 mm)

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

Eyepiece (Ø 23,2 mm): HWF 10 x / Ø 20 mm OBB-A1617

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OBB-A 16 17
Model-Serie	OBB-AE
Assortment	Optics
Product Type	Microscopes
Product group	Transmitted light microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761253759
Customs tariff number	90119000

FUNCTIONS

Standard



ENGLISH

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	85×65×45 mm
Gross weight	0.40 kg
Net weight	0.10 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	300 g
Delivery time	1 d