

KERN 348-08

Milligram weights, M1, 200 mg, Platelet, stainless steel polished



Article codes, product groups & Status information

Model article-code	348-08
Model-Serie	348
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761240209
Predecessor Model	348-8-2.022a
Planned start date	2020-1HJ
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary

Nominal value	200 mg
Tolerance (OIML)	0,6 mg
Case aluminium	317-009-600
Case plastic	347-009-400
Case wood	338-090-200

Technical data - Construction

Dimensions (W×D×H)	9×9×4 mm
Material	Stainless steel polished

Technical data - Measuring system

Nominal value	200 mg
OIML class	M1

Functions

Standard

OIML conform	Yes
Tolerance (OIML)	0,6 mg
Density (kg/m ³)	7.960

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	9×9×4 mm
Gross weight	0,05 kg
Net weight	0,05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkKS Calibration	962-658
-------------------	---------

KERN 348-08

Milligram weights, M1, 200 mg, Platelet, stainless steel polished



Option

