

PRODUCT BROCHURE

KERN 344-036

Date: 10.03.2026

ENGLISH



Weight set, M1, 1 g - 100 g, Knob, stainless steel fine turned, in aluminium case

CONTACT


KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	344-036
Model-Serie	344
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main article
EAN code (Model)	404576 1224872
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 g - 100 g
Case aluminium	314-030-600
Weight set case	in aluminium case

Technical data - Construction

Material	stainless steel fine turned
Is weight set	yes
Number of weights in set	9
Type of case	in aluminium case

Technical data - Measuring system

Nominal value	1 g - 100 g
OIML class	M1
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	165×125×75 mm
Gross weight	0.90 kg
Net weight	0.70 kg

Shipping method	Parcel service
Packaging component - by weight - plastics	200 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-616
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
962-616	DAkkS Calibration KERN 962-616
962-616R	DAkkS Recalibration KERN 962-616R
920-711	Marking (Last digits of serial number) KERN 920-711

SPARE PARTS

Model	Description
314-030-600	Aluminium Protective Box KERN 314-030-600
315-242	Tweezers KERN 315-242