

PRODUCT BROCHURE

KERN 314-064

Date: 09.03.2026

ENGLISH



Weight set, E2, 1 g - 1 kg, Knob, stainless steel polished, in plastic case

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	314-064
Model-Serie	314
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Main article
EAN code (Model)	4045761081062
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary

Nominal value	1 g - 1 kg
Case plastic	313-082-400
Weight set case	in plastic case

Technical data - Construction

Material	stainless steel polished
Is weight set	yes
Number of weights in set	13
Type of case	in plastic case

Technical data - Measuring system

Nominal value	1 g - 1 kg
OIML class	E2
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	350×320×180 mm
Gross weight	3.8 kg
Net weight	3.8 kg
Shipping method	Parcel service

Packaging component - 0 g
by weight - plastics

Delivery time 1 d

Technical data - Verification

Verification possible yes

Optional Services

DAkkS Calibration 962-319

FUNCTIONS

Standard



Option



SERVICES

Model	Description
962-319	DAkkS Calibration KERN 962-319
962-319R	DAkkS Recalibration KERN 962-319R
921-319	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-319
920-714	Marking (Last digits of serial number) KERN 920-714

SPARE PARTS

Model	Description
313-082-400	Plastic Weight Box KERN 313-082-400
315-242	Tweezers KERN 315-242