

PRODUCT BROCHURE

KERN 337-13

Date: 22.03.2025

ENGLISH



Individual weight, OIML Class F2, 5 kg, Knob, stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

- Individual weight, OIML Class F2, 5 kg, Knob, stainless steel fine turned (OIML)
- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	337-13
Model-Serie	337
Assortment	KERN
Product Type	Test weights
Product group	OIML F2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761006522
Customs tariff number	84239090

Technical data - Primary

Nominal value	5 kg
Tolerance (OIML)	80 mg
Case aluminium	317-130-600
Case plastic	347-130-400
Case wood	337-130-200

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	5 kg
OIML class	F2
OIML conform	yes
Tolerance (OIML)	80 mg
Density (kg/m ³)	7950
Volume (cm ³)	625

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	80×80×147 mm
Gross weight	5 kg
Net weight	5 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	no
-----------------------	----

Optional Services

DAkkS Calibration	962-443
-------------------	---------

FUNCTIONS

Standard



ACCESSORIES

Model	Description
317-280	Gloves KERN 317-280
317-290	Gloves KERN 317-290
337-130-200	Wooden Box for Set of Weights KERN 337-130-200
347-130-400	Plastic Weight Box KERN 347-130-400
317-130-600	Aluminium Protective Box KERN 317-130-600

SERVICES

Model	Description
962-443	DAkkS Calibration KERN 962-443
962-443R	DAkkS Recalibration KERN 962-443R
920-30	Marking with stroke number KERN 920-30
920-32	Marking by etching KERN 920-32
921-443	Guaranteed Plus-deviation ($\pm 1/3$ OIML) KERN 921-443
920-301	Marking with stroke number (Last digits of serial number) KERN 920-301
920-321	Marking by etching (last digits of the serial number) KERN 920-321

RELATED ARTICLES DATA

Model name	Nominal value	Tolerance (OIML)	Case aluminium	Case plastic	Case wood
337-01	1 g	0,3 mg	317-010-600	347-030-400	337-010-200
337-02	2 g	0,4 mg	317-020-600	347-030-400	337-020-200
337-03	5 g	0,5 mg	317-030-600	347-030-400	337-030-200
337-04	10 g	0,6 mg	317-040-600	347-050-400	337-040-200
337-05	20 g	0,8 mg	317-050-600	347-050-400	337-050-200
337-06	50 g	1 mg	317-060-600	347-070-400	337-060-200
337-07	100 g	1,6 mg	317-070-600	347-070-400	337-070-200
337-08	200 g	3 mg	317-080-600	347-080-400	337-080-200
337-09	500 g	8 mg	317-090-600	347-090-400	337-090-200

Model name	Nominal value	Tolerance (OIML)	Case aluminium	Case plastic	Case wood
337-11	1 kg	16 mg	317-110-600	347-110-400	337-110-200
337-12	2 kg	30 mg	317-120-600	347-120-400	337-120-200
337-13	5 kg	80 mg	317-130-600	347-130-400	337-130-200
337-14	10 kg	160 mg	317-140-600	347-140-400	337-140-200
337-141	10 kg	160 mg	not available	not available	337-141-200
337-15	20 kg	300 mg	317-150-600	-	337-150-200
337-151	20 kg	300 mg	not available	not available	337-151-200
337-16	50 kg	800 mg	317-160-600	-	337-160-200
337-161	50 kg	800 mg	not available	not available	337-161-200
337-131	5 kg	80 mg	-	-	337-131-200