


KERN 346-87

Individual weight, OIML Class M1, 10 kg, Block, cast iron lacquered



Article codes, product groups & Status information

Model article-code	346-87
Model-Serie	346
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761361799
Predecessor Model	346-87-2019a; 346-87-2024a
Successor Model	346-87-2024e2
Planned start date	2025-2HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	10 kg
Tolerance (OIML)	500 mg
Case aluminium	346-070-600

Technical data - Construction

Dimensions (W×D×H)	193×95×108 mm
Material	Cast iron lacquered

Technical data - Measuring system

Nominal value	10 kg
OIML class	M1
OIML conform	Yes
Tolerance (OIML)	500 mg

Functions

Standard

Density (kg/m ³)	7100
Volume (cm ³)	1408,4507

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	215×110×120 mm
Gross weight	10 kg
Net weight	10 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	Yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-644
-------------------	---------

KERN 346-87

Individual weight, OIML Class M1, 10 kg, Block, cast iron lacquered

