

PRODUCT BROCHURE

KERN 333-046

Date: 10.03.2026

ENGLISH



Weight set, F2, 1 mg - 200 g, Knob, stainless steel fine turned, in aluminium case

CONTACT


KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



DESCRIPTION

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	333-046
Model-Serie	333
Assortment	KERN
Product Type	Test weights
Product group	OIML F2 Test weights
Product usage type	Main article
EAN code (Model)	404576 1224544
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 mg - 200 g
---------------	--------------

Case aluminium 313-040-600

Weight set case in aluminium case

Technical data - Construction

Material stainless steel fine turned

Is weight set yes

Type of case in aluminium case

Technical data - Measuring system

Nominal value 1 mg - 200 g

OIML class F2

OIML conform yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	190×130×80 mm
Gross weight	1.2 kg
Net weight	0.90 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	200 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-403
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
962-403	DAkkS Calibration KERN 962-403
962-403R	DAkkS Recalibration KERN 962-403R
921-403	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-403
920-704	Marking (Last digits of serial number) KERN 920-704

SPARE PARTS

Model	Description
313-040-600	Aluminium Protective Box KERN 313-040-600