

KERN 357-04

Individual weight, OIML Class M2, 10 g, Knob, stainless steel fine turned



Article codes, product groups & Status information

Model article-code	357-04
Model-Serie	357-04
Assortment	KERN
Product Type	Test weights
Product group	OIML M2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761288959
Planned start date	2022-1HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	10 g
Tolerance (OIML)	6 mg
Case aluminium	317-040-600
Case plastic	347-050-400
Case wood	337-040-200

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	10 g
OIML class	M2
OIML conform	Yes
Tolerance (OIML)	6 mg

Functions

Standard

Density (kg/m ³)	7.950
Volume (cm ³)	1,1905

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	24×24×51 mm
Gross weight	0,05 kg
Net weight	0,05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkKS Calibration	962-634
-------------------	---------

KERN 357-04

Individual weight, OIML Class M2, 10 g, Knob, stainless steel fine turned



Option

