

## PRODUCT BROCHURE

---

# OPTICS OBB-A1210

Date: 30.04.2025

ENGLISH



Trinocular lens barrel

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## DESCRIPTION

Microscope tube OBB-A1210

## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                                 |
|-----------------------|---------------------------------|
| Model article-code    | OBB-A1210                       |
| Model-Serie           | OBB-AT                          |
| Assortment            | Optics                          |
| Product Type          | Microscopes                     |
| Product group         | Polarising microscopes          |
| Product usage type    | Spare part;Accessories/Services |
| EAN code (Model)      | 4045761167421                   |
| Customs tariff number | 90119000                        |

### Technical data - Primary

|      |            |
|------|------------|
| Tube | Trinocular |
|------|------------|

### Technical data - Construction

|                    |                |
|--------------------|----------------|
| Tube type          | Siedentopf     |
| Tube angle         | 30°            |
| Tube rotation      | yes            |
| Eye distance (max) | 75 mm          |
| Eye distance (min) | 50 mm          |
| Adjustable dioptré | one-sided      |
| Dimensions (W×D×H) | 230×120×100 mm |

### Technical data - Environmental conditions

|                  |       |
|------------------|-------|
| Storage Temp Min | -5 °C |
| Stotage Temp Max | 40 °C |

### Technical data - Packing & Shipping

|                                            |                |
|--------------------------------------------|----------------|
| Dimensions packaging (W×D×H)               | 265×185×165 mm |
| Gross weight                               | 2 kg           |
| Net weight                                 | 1.4 kg         |
| Shipping method                            | Parcel service |
| Packaging component - by weight - plastics | 600 g          |
| Delivery time                              | 1 d            |

## FUNCTIONS

|          |
|----------|
| Standard |
|----------|

