

KERN OZB-A4571UK

Professional illumination units guarantee outstanding, uniform and strong illumination



Article codes, product groups & Status information

Model article-code	OZB-A4571UK
Model-Serie	OZB-IR
Assortment	Optics
Product Type	Microscopes
Product group	External Illumination
Product usage type	Main article
EAN code (Model)	4045761218482
Planned start date	2019-1HJ
Customs tariff number	94054239

Technical data - Primary

Illumination type	LED
Illumination intensity	4W

Technical data - Construction

Color temperature	7000-11000 K
Dimmable	Yes
Ring diameter	6 cm

Technical data - Power supply

Protection Class	Class III - Low voltage
------------------	-------------------------

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

Technical data - Packing & Shipping

Functions

Standard

Dimensions packaging (W×D×H)	235×200×50 mm
------------------------------	---------------

Gross weight	0.50 kg
--------------	---------

Net weight	0.45 kg
------------	---------

Shipping method	Parcel service
-----------------	----------------

Packaging component - by weight - paper	0 g
---	-----

Packaging component - by weight - cardboard	72.0324 g
---	-----------

Packaging component - by weight - plastics	2.3112 g
--	----------

Packaging component - by weight - styrofoam	21.9564 g
---	-----------

Packaging component - by weight - foam	0 g
--	-----

Delivery time	1 d
---------------	-----

Technical data - Verification

WEEE	Yes
------	-----

KERN OZB-A4571UK

Professional illumination units guarantee outstanding, uniform and strong illumination

