

PRODUCT BROCHURE

KERN 317-05-2022i

Date: 15.03.2025

ENGLISH



Individual weight, OIML Class E2, 20 g, Knob, stainless steel polished

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	317-05-2022i
Model-Serie	317
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Trade parts list element
EAN code (Model)	4045761245310
Predecessor Model	317-05; 317-05-2020a
Successor Model	317-05
Customs tariff number	84239090

Technical data - Primary

Nominal value	20 g
Tolerance (OIML)	0,08 mg
Case aluminium	317-050-600
Case plastic	317-050-400
Case wood	317-050-100

Technical data - Construction

Material	stainless steel polished
----------	--------------------------

Technical data - Measuring system

OIML class	E2
OIML conform	yes
Density (kg/m ³)	7950
Volume (cm ³)	2,5

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	13×13×22 mm
Gross weight	0.05 kg
Net weight	0.05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-335
-------------------	---------

FUNCTIONS

Standard



ACCESSORIES

Model	Description
317-050-100	Wooden Box for Set of Weights KERN 317-050-100
317-050-400	Plastic Weight Box KERN 317-050-400
317-050-600	Aluminium Protective Box KERN 317-050-600
315-242	Tweezers KERN 315-242

SERVICES

Model	Description
952-335	Verification KERN 952-335
962-335	DAkkS Calibration KERN 962-335
962-335R	DAkkS Recalibration KERN 962-335R
920-33	Marking by etching KERN 920-33
921-335	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-335
920-331	Marking by etching (last digits of the serial number) KERN 920-331