

PRODUCT BROCHURE

KERN 313-066

Date: 10.03.2026

ENGLISH



Weight set, E2, 1 mg - 1 kg, Knob, stainless steel polished, in aluminium case

CONTACT


KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



DESCRIPTION

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	313-066
Model-Serie	313
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Main article
EAN code (Model)	4045761223547
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary

Nominal value	1 mg - 1 kg
---------------	-------------

Case aluminium	313-060-600
Weight set case	in aluminium case

Technical data - Construction

Material	stainless steel polished
Is weight set	yes
Type of case	in aluminium case

Technical data - Measuring system

Nominal value	1 mg - 1 kg
OIML class	E2
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	230×140×110 mm
Gross weight	3.6 kg
Net weight	2.8 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	800 g
Delivery time	1 d

Technical data - Verification

Verification possible	yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-305
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
962-305	DAkkS Calibration KERN 962-305
962-305R	DAkkS Recalibration KERN 962-305R
921-305	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-305
920-706	Marking (Last digits of serial number) KERN 920-706

SPARE PARTS

Model	Description
<u>313-060-600</u>	Aluminium Protective Box KERN 313-060-600